

全基因组的比较基因组杂交技术介绍

(Whole-Genome and Custom Fine-Tiling Array CGH)

Whole-Genome and Custom Fine-Tiling Array CGH

Comparative Genomic Hybridization (CGH) measures DNA copy number differences between a reference genome and your sample genome. NimbleGen offers two high-definition array CGH products: whole-genome array CGH and fine-tiling array CGH. Whole-genome CGH measures copy number differences in DNA across entire genomes, while fine-tiling CGH determines breakpoints at ultrahigh resolution in DNA targets located in any subset of DNA (either a contiguous region or many different regions of interest).

Key features of NimbleGen's microarrays - ultra-high density and long oligo probes - yield a powerful, high-definition CGH platform for whole genome analysis. With the flexibility of NimbleGen array synthesis, chromosomal aberrations identified in whole-genome surveys can be readily examined at even higher resolution with custom fine-tiling array CGH. Such studies can map genomic changes down to the gene and even exon level and enable breakpoint mapping to less than 500 bp intervals. NimbleGen's array CGH detects both single and multiple copy changes.

Advantages of NimbleGen Array CGH

Highest Resolution Array CGH

NimbleGen's high capacity microarrays contain 385,000 probes on a single glass slide. This feature density enables design of a human whole genome CGH array with probes tiled through genic and intergenic regions at a median probe spacing of 6,000 bp. This tiling path array design approach ensures you get the complete picture on genome-wide copy number changes. Ultra-high resolution detection of small deletions and amplifications, as well as breakpoint mapping to less than 500 bp intervals, can be achieved using NimbleGen's fine-tiling arrays with probes as dense as 10 bp spacing.

Isothermal (T_m-Balanced), Long-Oligonucleotide Probes

NimbleGen's exclusive isothermal array design approach enables uniform probe performance, eliminating hybridization artifacts and/or bias and providing higher quality data. Probe lengths are adjusted (45mer - 85mer) to equalize the melting

temperature ($T_m = 76^\circ\text{C}$) across the entire set. Thus, probes are optimized to perform equivalently at a given stringency in the genomic regions your experiments require, including AT- and GC-rich regions.

Custom Array CGH

NimbleGen's highly flexible Maskless Array Synthesis technology enables rapid prototyping of new array designs. This flexibility enables cost-effective array CGH design for your organism of interest at the whole-genome level or focused on specific chromosomal regions.

Most Up-to-date Genome Build

NimbleGen CGH arrays are always current with the latest genome builds, so there's no need to put your experiments on hold for a new array CGH design. In addition, you can continue to access array designs based on past genome data builds, which can be particularly useful for comparisons to prior studies.

Products and Services

Human Whole-Genome Array CGH

With a capacity of 385,000 isothermal probes, NimbleGen's high-density, long oligo arrays span the entire non-repetitive regions of the human genome in a single array, tiling the full genome at a median probe spacing of 6,000 bp. NimbleGen's Human Whole-Genome Array CGH, available as a catalog design (no design fees), is a cost-effective platform for genome-wide analysis of copy number changes. Unlike other commercial whole-genome array CGH platforms, NimbleGen's Human Whole-Genome Array CGH tiles through both genic and intergenic regions of the genome, providing the most thorough, unbiased coverage available.

Human Whole-Genome Array CGH Analysis of Human Tumor Cell Line

The large array capacity (385,000 features), long oligo probes, and tiling path microarray design enables genome-wide detection of copy number changes in genic and intergenic regions.

Figure A: Whole-genome array CGH results displayed in an “all chromosomes” view in NimbleGen's SignalMap browser, which includes gene annotation with direct links to NCBI. Data are provided as General Feature Format (GFF) files for use with

SignalMap and other genome browsers. [[Click for SignalMap view](#)]





Figure B: Whole-genome plot of data processed with a 50 Kb averaging window for a global view of copy number changes.

Figure C: Zoomed-in view of a region in chromosome 3 showing a single copy loss and 50 Kb microdeletion. Top track is gene annotation, middle track is averaged data with segmentation analysis (red horizontal lines demarcate copy number changes), and bottom track is individual probe data (unaveraged).

Custom Fine-Tiling Array CGH

NimbleGen's Custom Fine-Tiling Array CGH can be used to detect deletions and amplifications and to map the associated breakpoints with unprecedented resolution. You can choose your region(s) of interest for a fine-tiling array design with probes spaced as dense as 10 bp. This density of probe placement enables ultra-high resolution mapping of breakpoints to an interval that can be validated by PCR amplification and sequencing. The selected regions need not be contiguous or even within the same chromosome. For example, the copy number changes identified in a whole-genome array CGH experiment can be combined and tiled at higher resolution on a Custom Fine-Tiling CGH Array.



Fine-Tiling Array CGH Analysis of Human Tumor Cell Line

Regions from four separate chromosomes were tiled at high resolution, enabling precise determination of breakpoints.

Figure D: Full fine-tiling data set in SignalMap browser

Figure E: Detection of an ~80 Kb deletion

Figure F: Detection of a single-copy change breakpoint

Applications

NimbleGen array CGH is a significant advance over current array CGH methods for the detection of DNA copy number changes and mapping of the associated breakpoints. While these changes are commonly recognized as the underlying basis for many congenital disorders and complex diseases (such as cancer), variation in copy number is also commonly found between healthy individuals. Characterizing normal versus aberrant chromosomal differences will be key to understanding the molecular mechanisms for a given disease.

NimbleGen array CGH provides the highest level resolution for detecting a wide range of copy number changes on either a genome-wide or fine-tiling level, including:

Homozygous and hemizygous deletions

Single and multiple copy amplifications

Unbalanced translocations

NimbleGen Service Overview

NimbleGen's CGH microarray service consists of the following steps:

Customer and NimbleGen create custom array design, and NimbleGen manufactures the array.

Customer supplies 1-3 μ g purified, undegraded sample and reference genomic DNA (reference optional).

NimbleGen labels the sample, performs the CGH hybridization, scans the array, extracts the data, and performs the segmentation analysis.

Data Delivery

The data delivered with NimbleGen's CGH service includes:

Raw data

Segmentation analysis based on the DNACopy package

Segmentation GFF file

Genome annotation GFF file

SignalMap™ GFF visualization software

Data sets are provided in GFF file format for easy data set analysis with the SignalMap™ data browser software. SignalMap displays the segmentation analysis in relation to the known state of sequence annotation, including genes, exons, and repetitive elements.

NimbleGen Array CGH General Specifications

CGH Array Specifications

Probe length	Isothermal (T _m =76°C, 45mer - 85mer)
Probe design format	Tiled throughout genome on forward strand. Repeat sequences masked.
Total features	385,000
Feature Size	16 μ m x 16 μ m
Array size	17.4mm x 13mm
Slide size	1" x 3" (25mm x 75mm) glass
Sample required / array	1 - 3 μ g isolated gDNA test sample 1 - 3 μ g isolated gDNA reference sample

Sample labeling	2 color (Cy3 and Cy5)
Sequence source	UCSC Genome Browser - http://genome.ucsc.edu/
Deliverables	Raw and processed data, segmentation analysis, SignalMap™ data browser

Available Array Designs

Whole Genome CGH Arrays

Available	Human Tiled	Mouse Tiled	C. elegans Tiled
Genomes	Whole Genome	Whole Genome	Whole Genome
Genome Build	HG17	MM5	ce2, WS120
Median Probe Spacing	6000 bp	5000 bp	140 bp
Probe Length	isothermal probes: Tm 76 C, length 45-85 bp.		

Fine-Tiling CGH Arrays

Available	Human CHR X	HG17
Genomes	Fine-Tiling	
Median Probe Spacing	106 bp	
Probe Length	isothermal probes: Tm 76 C, length 45-85 bp	

Custom Fine-Tiling CGH Arrays

Available	Any genome for which high-quality sequence is available.
Genomes	
Median Probe Spacing	Researcher Specified from 10bp to 5000bp
Bioinformatics	Custom array design available for an additional fee

Resolution	Custom tiling up to 385,000 probes per array. Single or multiple chromosomal regions per array. Single or multiple array designs.
------------	---
