

Caustic 苛性碱

订货号	包装	英文名称	中文名称
1.05012.1000	1kg	Potassium hydroxide pellets pure	超纯级氢氧化钾
1.05012.5000	5kg	Potassium hydroxide pellets pure	超纯级氢氧化钾
1.05012.9050	50kg	Potassium hydroxide pellets pure	超纯级氢氧化钾
1.05033.0500	500g	Potassium hydroxide pellets GR for analysis	优级纯氢氧化钾
1.05033.1000	1kg	Potassium hydroxide pellets GR for analysis	优级纯氢氧化钾
1.05033.5000	5kg	Potassium hydroxide pellets GR for analysis	优级纯氢氧化钾
1.05033.9025	25kg	Potassium hydroxide pellets GR for analysis	优级纯氢氧化钾
1.05033.9050	50kg	Potassium hydroxide pellets GR for analysis	优级纯氢氧化钾
1.05029.1000	1kg	Potassium hydroxide pellets GR for analysis (max. 0,05% Na) ACS,Reag. Ph Eur	中钠氢氧化钾
1.05029.9050	50kg	Potassium hydroxide pellets GR for analysis (max. 0,05% Na) ACS,Reag. Ph Eur	中钠氢氧化钾
1.05021.0250	250g	Potassium hydroxide pellets GR for analysis (max. 0,002% Na) ACS,ISO,Reag. Ph Eur	低钠氢氧化钾
1.05021.1000	1kg	Potassium hydroxide pellets GR for analysis (max. 0,002% Na) ACS,ISO,Reag. Ph Eur	低钠氢氧化钾
1.05021.5000	5kg	Potassium hydroxide pellets GR for analysis (max. 0,002% Na) ACS,ISO,Reag. Ph Eur	低钠氢氧化钾
1.05021.9025	25kg	Potassium hydroxide pellets GR for analysis (max. 0,002% Na) ACS,ISO,Reag. Ph Eur	低钠氢氧化钾
1.05545.1000	1L	Potassium hydroxide solution 47% GR for analysis	优级纯氢氧化钾溶液
4.80311.1000	1L	Potassium hydroxide solution 40% (max. 500 ppm Na) for analysis	低钠氢氧化钾溶液
4.80311.9025	25L	Potassium hydroxide solution 40% (max. 500 ppm Na) for analysis	低钠氢氧化钾溶液
1.06462.1000	1kg	Sodium hydroxide pellets pure	超纯级氢氧化钠
1.06462.5000	5kg	Sodium hydroxide pellets pure	超纯级氢氧化钠
1.06462.9050	50kg	Sodium hydroxide pellets pure	超纯级氢氧化钠
1.06498.0500	500g	Sodium hydroxide pellets GR for analysis ISO	ISO级氢氧化钠
1.06498.1000	1kg	Sodium hydroxide pellets GR for analysis ISO	ISO级氢氧化钠
1.06498.5000	5kg	Sodium hydroxide pellets GR for analysis ISO	ISO级氢氧化钠
1.06498.9025	25kg	Sodium hydroxide pellets GR for analysis ISO	ISO级氢氧化钠
1.06498.9050	50kg	Sodium hydroxide pellets GR for analysis ISO	ISO级氢氧化钠
1.06469.1000	1kg	Sodium hydroxide pellets GR for analysis (max. 0,02% K) ACS,Reag. Ph Eur	中钾氢氧化钠
1.06469.5000	5kg	Sodium hydroxide pellets GR for analysis (max. 0,02% K) ACS,Reag. Ph Eur	中钾氢氧化钠
1.06469.9050	50kg	Sodium hydroxide pellets GR for analysis (max. 0,02% K) ACS,Reag. Ph Eur	中钾氢氧化钠
1.06495.0250	250g	Sodium hydroxide pellets GR for analysis (max. 0,0002% K) ACS,Reag. Ph Eur	低钾氢氧化钠
1.06495.1000	1kg	Sodium hydroxide pellets GR for analysis (max. 0,0002% K) ACS,Reag. Ph Eur	低钾氢氧化钠
1.06495.5000	5kg	Sodium hydroxide pellets GR for analysis (max. 0,0002% K) ACS,Reag. Ph Eur	低钾氢氧化钠
1.06495.9050	50KG	Sodium hydroxide pellets GR for analysis (max. 0,0002% K) ACS,Reag. Ph Eur	低钾氢氧化钠
1.06466.0050	50g	Sodium hydroxide monohydrate 99,99 Suprapur	单水合氢氧化钠
1.06466.0500	500g	Sodium hydroxide monohydrate 99,99 Suprapur	单水合氢氧化钠
1.06467.9010	10kg	Sodium hydroxide granulated	粒状氢氧化钠
1.06467.9050	50kg	Sodium hydroxide granulated	粒状氢氧化钠
1.05587.2500	2.5L	Sodium hydroxide solution about 32% extra pure	超纯级氢氧化钠溶液
1.05587.9025	25L	Sodium hydroxide solution about 32% extra pure	超纯级氢氧化钠溶液
1.05587.9200	200L	Sodium hydroxide solution about 32% extra pure	超纯级氢氧化钠溶液
1.05588.1000	1L	Sodium hydroxide solution min. 10% (1.11) GR for analysis	10%氢氧化钠溶液
1.05588.9010	10L	Sodium hydroxide solution min. 10% (1.11) GR for analysis	10%氢氧化钠溶液
1.05589.0250	250ml	Sodium hydroxide solution 30% Suprapur	32%氢氧化钠溶液
1.05589.1000	1L	Sodium hydroxide solution 30% Suprapur	32%氢氧化钠溶液
1.05590.2500	2.5L	Sodium hydroxide solution about 32% for analysis (for the determination of nitrogen)	氮滴定氢氧化钠溶液
1.05590.9025	25L	Sodium hydroxide solution about 32% for analysis (for the determination of nitrogen)	氮滴定氢氧化钠溶液
1.05591.2500	2.5L	Sodium hydroxide solution min. 27% (1.30) GR for analysis (for the determination of nitrogen)	27%氢氧化钠溶液
1.05591.9025	25L	Sodium hydroxide solution min. 27% (1.30) GR for analysis (for the determination of nitrogen)	27%氢氧化钠溶液
1.05593.9025	25L	Sodium hydroxide solution 21% GR for analysis	21%氢氧化钠溶液
1.11360.2500	2.5L	Sodium hydroxide solution min. 45% GR for analysis	45%氢氧化钠溶液
1.11360.9025	25L	Sodium hydroxide solution min. 45% GR for analysis	45%氢氧化钠溶液
1.58793.1000	1L	Sodium hydroxide solution 50% GR for analysis	50%氢氧化钠溶液
1.58793.9025	25L	Sodium hydroxide solution 50% GR for analysis	50%氢氧化钠溶液